

## CENTRIFUGES

### Medico centrifuges

Rotek Plus offers a modified model of medico centrifuge which is ideal for routine centrifuging tests of medical practitioners, hospitals, colleges, institutions and industrial laboratories etc. Fitted with resiliently mounted motor for vibration free performance. A specially wound 5 step transformer is used for regulation. The highly sophisticated elegant body of the equipment is made with fibre reinforced plastic which gives protection against corrosion and electric shock.

Outer body made of aluminium and powder coated model also available.

Maximum speed 4000RPM

Maximum RCF 1880 x g

Motor 1/20 HP AC/DC universal

Specifications	Cat No.	Capacity
MEDICO Centrifuge - Swing out head(Outer body made of fibre glass and paint finished)	PMC - 4F	4 x 15 ml
	PMC - 6F	6 x 15 ml
Medico Centrifuge - Swing out head	PMC - 4 A	4 x 15 ml
(Outer body made of aluminium and powder coated)	PMC - 6 A	6 x 15 ml



### CLINICAL CENTRIFUGE

Clinical Centrifuge is specially designed for the centrifuging requirements of medical practitioners and pathological laboratories, hospitals and such other institutions. The fixed angle rotor head is made of aluminium alloy, fitted with resiliently mounted motor for vibration free performance. A specially wound 5 step transformer speed regulator is used for speed control.



CLINICAL	PCH - 4	4 x 15 ml
Clinical Centrifuge - Angle rotor head	PCH -6	6 x 15 ml
	PCH -8	8 x 15 ml

## MULTI PURPOSE LABORATORY CENTRIFUGE

Multipurpose laboratory centrifuges are characterized by extreme reliability and a long working life, ideal for medical, institutional and industrial laboratories. Special feature includes stepless electronic speed regulator, sturdy body etc.

specifications	Cat.No.
Multipurpose Laboratory Centrifuge with speed regulator(without head)	P 8 A
Multipurpose Laboratory Centrifuge with 30 min spring Timer and speed regulator (without head)	P 8 B
Above with Digital Timer	P8 BE
Multipurpose Laboratory Centrifuge with Digital Timer and digital RPM indicator (without head)	P8 E



## ROTOR HEADS FOR P 8 A, P 8 B, P 8 E, P8BE (INTERCHANGABLE)

	Cat. No.	Capacity in ml.
Swing out heads	SH - 1	8 x 15
	SH - 2	4 x 50
Angle rotor heads	AH - 1	8 x 15
	AH - 2	12 x 15
	AH - 3	4 x 50
	AH - 4	6 x 50
	AH - 5	16 x 15