

DEEP FREEZER (HORIZONTAL / VERTICAL) – 20°C

Double walled construction outer body made of G.I. sheet, powder coated, inner chamber made of stainless steel, insulation PUF . The outer door is provided with magnetic , beading, lock, key and handle, temp: range -20oC, fitted with digital temp: controller cum indicator with built in on delay . Equipped with CFC free compressor and refrigerant. The entire unit is mounted on castor wheels and lock. Vertical models are supplied with stainless steel perforated trays.

Outer body with stainless steel 304 quality sheet with GMP standard also available.

Non use of stabilizer will render the warranty of the equipment invalid.

Stabilizer rating:- 1 KVA/ 2KVA / 3KVA.

HORIZONTAL

Cat. No.	Size cms			Capacity
	L	W	H	
PDF - 110	60	45	45	110 ltr
PDF - 175	86	51	51	175 ltr
PDF - 280	120	53	53	280 ltr
PDF - 450	155	55	55	450 ltr.

VERTICAL

Cat. No.	Size cms			Capacity
	L	W	H	
PDFV - 100	45	45	60	100 ltr
PDFV - 175	51	51	67	175 ltr
PDFV - 280	53	53	90	280 ltr
PDFV - 400	65	63	98	400 ltr.





DEEP FREEZER (HORIZONTAL / VERTICAL)-35°C

Double walled construction, outer body made of G.I. sheet, powder coated, inner chamber made of stainless steel, insulation PUF. The outer door is provided with magnetic beading, lock, key and handle, temp: 35 o C fitted with digital temp: controller cum indicator with built in on delay. Equipped with CFC free compressor and refrigerant.

The entire unit is mounted on castor wheels and lock

Vertical models are supplied with stainless steel perforated trays.

Outer body with stainless steel 304 quality sheet with GMP standard also available

Suitable voltage stabilizer is to be used to protect the equipment from voltage fluctuations.

Non use of stabilizer will render the warranty of the equipment invalid.

Stabilizer rating :- 1 KVA /2 KVA /3 KVA

HORIZONTAL

Cat. No.	Size cms	Capacity
	L W H	
PDFL - 110	60x45x45	110 ltr
PDFL - 175	86x51x51	175 ltr
PDFL - 280	120x53x53	280 ltr

VERTICAL

Cat. No.	Size cms	Capacity
	L W H	
PDFVL - 100	45x45x60	100 ltr
PDFVL- 175	51x51x67	175 ltr
PDFVL- 280	53x53x90	280 ltr





DEEP FREEZER (HORIZONTAL / VERTICAL) - 80°C (Ultra Low Temp.)

Double walled construction, outer body made of G.I. sheet, inner chamber made of stainless steel, insulation PUF. The outer door is provided with magnetic beading, lock, key and handle, temperature -80°C (at ambient temp: of 20°C AC room) fitted with digital temp: controller cum indicator with special features of cascade function and built in on delay for compressor. Equipped with CFC free IMPORTED compressor and refrigerant. The entire unit is mounted on castor wheels and lock. Vertical models are supplied with stainless steel perforated trays.

Outer body with stainless steel 304 quality sheet with GMP standard also available.

Suitable voltage stabilizer is to be used to protect the equipment from voltage fluctuations.

Non use of stabilizer will render the warranty of the equipment invalid.

Stabilizer rating :- 4 K.V.A / 5 K.V.A (Servo Stabilizer)

HORIZONTAL

Cat. No.	Size cms			Capacity
	L	W	H	
PUDF - 110	60	45	45	110 ltr
PUDF - 175	67	51	51	175 ltr
PUDF- 250	90	53	53	250 ltr

VERTICAL

Cat. No.	Size cms			Capacity
	L	W	H	
PUDFV - 100	45	45	60	100 ltr
PUDFV- 175	51	51	67	175 ltr
PUDFV- 280	53	53	90	280 ltr

