



DEPTH SAMPLER

Designed to collect water samples from desired depth, suitable to collect water from pond, lake, and sea etc, easy to operate. The sample collection container is made of transparent Acrylic tube. Provided with drain cock in the bottom, provided with shutter opening and lid locking arrangement, one messenger is made of M.S. and galvanized is supplied after lowering the depth sampler to the desired depth, the messenger can be send down through the rope, that will hit on the lock and release the lock simultaneously both the shutters will close and water in that depth will be stored in the container. Supplied with 15 mtr rope.

Cat No.	Capacity
PDS - 1	2 ltr.