

DRY BLOCK DIGESTOR

These dry block digestors are suitable for Digestion, Boiling, and various industrial applications, heating block is made with cast aluminium alloy, suitable to work up to 400°C. Block provided with 24 nos holes with 28 mm dia and 120mm depth suitable to insert 25 mm dia 200 mm long digestion tube. Heaters are placed in the base of the block. Block mounted on heavy M.S. sheet frame, properly painted. Model available with thermostat upto 300°C and digital temp: controller model for 400°C.



Cat No.	Temp Range
PDBD -01	50°C to 300°C
PDBD -02	50°C to 400°C

Digital Temp. controller with Soak timer can be supplied on extra cost