

HEATING MANTLES

Suitable for laboratories of chemicals, pharmaceuticals, dyes, glass plants and distillation units etc., particularly useful for heating inflammable liquids. Mantle is made from fibre glass yarn and good quality nichrome/kanthal wire element. The fibre glass yarn can withstand temp: up to 350°C, mounted in a powder coated Aluminium body with glass wool insulation. Fitted with Energy regulator. Suitable to work on 220 VAC supply.



Cat No.	Capacity of flasks	Watts
PHM - 1	250 ml	120
PHM - 2	500 ml	200
PHM - 3	1 ltr.	300
PHM - 4	2 ltr.	500
PHM - 5	5 ltr.	800
PHM - 6	10 ltr.	1000
PHM - 7	20 ltr.	1500

SOXHELT EXTRACTION MANTLE

Ideal apparatus for the soxhelt and other types of extraction purpose in laboratory . Heating mantles arranged in a row. Housing is made with M.S. Sheet, duly powder coated, supplied with stand and adjustable rod. Supplied without any glass parts. Suitable to work on 220 VAC supply.



Cat No	Capacity	No. of flasks
PHMS-3S	250 ml	3
PHMS-4S	250 ml	4
PHMS- 6S	250 ml	6
PHMS - 3M	500 ml	3
PHMS - 4M	500 ml	4
PHMS - 6M	500 ml	6
PHMS - 3L	1000 ml	3
PHMS - 4L	1000 ml	4
PHMS - 6L	1000 ml	6