

HOT PLATE (ROUND)

Aluminium / Cast iron hot plate, available with different dia and different rating mounted on powder coated M.S. sheet fabricated body, available with energy regulator and without energy regulator. Maximum surface temp: 400°C

Cat. No.	Hot plate / dia in cms	
PHP - 1C	Cast iron Top	20
PHP - 1 A	Alu. Top	20
PHP - 1S	S. S Top	20
PHP - 2 A	Alu. Top	25
PHP - 2S	S. S Top	25



Thermostat 50°C to 300°C (imported) can be provided on extra cost (Accuracy +/- 5° C) Digital Temp: controller can be fitted at an extra cost (Accuracy +/- 5° C)

HOT PLATE (RECTANGULAR)

M.S. powder coated body. Hot plate surface is fabricated with thick SS sheet. Temp: regulated by energy regulator, available in different sizes and different rating.

Maximum surface temp: 400°C

Cat. No.	Size of Hot Plate	Rating
PHPR -1	25 x 30 cms.	1.2 KW
PHPR -2	25 x 40 cms.	1.6 KW
PHPR -3	30 x 45 cms.	2.0 KW

Thermostat 50°C to 300°C (imported) can be provided on extra cost

