

INSPIRATOR

This is an ideal equipment for preparing the media, inner chamber made of copper sheet of high thermal conductivity. The double walled water jacket is made of copper sheet, housed in wooden body, insulated with glass wool, the inner surface and tube holders are black epoxy coated for even radiation and distribution temp: in the chamber. Provided an acrylic framed top. Provision for water inlet and out let, the out let is the over flow tube. Both the inlet and over flow nipple are kept at same level so either one can be used as inlet or outlet, a drain valve is provided at the bottom .

Theater is provided at the outside bottom of the jacket, temperature is controlled by imported thermostat. Suitable to work on 220 VAC supply.

The chamber can accommodate 100 tube of media at time.



Cat. No.	Heater
PIN -1	500 watts