

KHAN SHAKER (Reciprocating Type)

Reciprocating action, forward and backward movement, available in two models, pre determined speed shaker (150 strokes per minute) and controllable speed shaker with AC TEFC motor. Tray with lotus clamp or rubber disc supporters can accommodate Erlenmeyer flasks and Tray with adjustable roller supporters can accommodate Khan racks, bottle of different sizes, displacement approx 40 mm. Various sizes are available.



Cat No.	Platform Capacity	Platform Size in cm
PKS - 01	9 flasks of 250 ml	38 x 38
PKS - 02	12 flasks of 250 ml	38 x 50
PKS - 03	16 flasks of 250 ml	50 x 50

While placing the order please specify the platform requirement.

KHAN SHAKER WITH CONTROLLABLE SPEED

Cat No.	Platform Capacity	Platform Size in cm
PKS - 04	9 flasks of 250 ml	38 x 38
PKS - 05	12 flasks of 250 ml	38 x 50
PKS - 06	16 flasks of 250 ml	50 x 50

While placing the order please specify the platform requirement