

MUFFLE FURNACE

Light weight construction, insulated with ceramic wool to reduce heat loss and weight. Maximum temperature 1000°C working temp: 900 °C, controlled by Energy Regulator, fitted with Thermo - couple, pyrometer Analog type. Provided with thermal fuse on specific request.

Model with Digital Display temp. controller cum indicator is also available.

Cat No.	Internal dimension in cm	Kilo watts
PMF - 1	10 x 10 x 22	1.5
PMF - 2	12 x 12 x 25	2.0
PMF - 3	15 x 15 x 22	2.5
PMF - 4	15 x 15 x 30	3.0



HIGH TEMP: MUFFLE FURNACE

High temperature muffle furnace with maximum temperature 1150°C, working temperature 1100°C. Smaller size muffle concealed type and bigger size muffle fabricated with grooved refractory, open element type.

Cat No.	Internal dimension in cm	Kilo watts
PMFH - 1	10 x 10 x 22	1.5
PMFH - 2	12 x 12 x 25	2.0
PMFH - 3	15 x 15 x 22	3.0
PMFH - 4	15 x 15 x 30	3.5