

SLIDE WARMING TABLE (THERMOSTATIC)

Thermostatic controlled slide warming table to maintain the temperature below the melting point of wax. Double walled construction, table top model, top plate made of choice of aluminium/stainless steel / Tin plated copper sheet mounted on a insulated outer body, fitted with imported thermostat, Suitable to work on 220 VAC supply.

Cat No.	Dimensions	Rating	Description
PSWT - 01	45 x 30 x 15 cm	(500 watts)	Top Al sheet
			Top S.S. sheet
			Top Copper sheet
PSWT - 02	60 x 45 x 15 cm	(1000 watts)	Top Al sheet
			Top S.S. sheet
			Top Copper sheet

