



**Raték Plus**

## TISSUE CULTURE RACK

Made with tubular M.S. pipe or tubular Aluminium pipe, duly powder coated. Having rack size of 120 cms length and 53 cms width, topped with thick plexi glass, height between each rack is 38 cms. Each rack is fitted with two fluo: tube and one U.V. Germicidal tube. Copper chokes / Electronic chokes are fitted on a separate box to avoid heat generation to the culture room (by keeping the choke box outside the room). The choke box is provided with cooling fan and castor wheel.



Cat. No.	No. of Racks	Material of Construction
PTCR - 01	4	MS Pipe & MS Frame
PTCR - 02	4	Aluminium pipe & Alu. Frame
PTCR - 03	5	MS Pipe & MS Frame
PTCR - 04	5	Aluminium pipe & Alu. Frame

Same with MS angle frame and topped with 6mm glass. 2 flourscent tubes provided in each rack.



AN ISO 9001-2015 CERTIFIED CO., REG: NO.02032018- AMAR-2018