

VACUUM OVEN RECTANGULAR

Double walled construction, outer body made with M.S sheet duly powder coated or painted. inner chamber made of stainless steel 304 grade sheet of thick guage to with stand vacuum, provided with ceramic wool and glass wool insulation to avoid thermal loss. Temperature range 50°C. to 150° C. Temperature controlled by imported hydraulic thermostat accuracy +/- 2°C. Digital temp: controller fitted model also available. Provided with valve and vacuum guage. Suitable to works on 220V AC Supply. Maximum vacuum 760 mm of mercury.

| Cat. No. | Inner Chamber Dimension |
|----------|-------------------------|
| PVO-R -2 | 35 X 35 X 35 |
| PVO-R-3 | 45 X 45 X 45 |
| PVO-R-4 | 45 X 45 X 60 |

Supplied without vacuum pump.

VACUUM OVEN ROUND

Double walled construction, outer body made with M.S. sheet duly powder coated / painted. Inner chamber made of stainless steel, glass wool insulation, temperature range 50°C. to 150°C. Temperature controlled by imported hydraulic thermostat / digital temp: controller. Provided with valve and vacuum gauge, see through door with toughened glass. Suitable to work on 230 V AC. Supplied without vacuum pump. Maximum vacuum 760 mm of mercury. Supplied without vacuum pump.



| Cat. No. | Chamber Size | Watts |
|----------|----------------------------|-------|
| PVO - 9 | 22.5 cms dia x 30 cm depth | 600 |
| PVO - 12 | 30 cms dia x 38 cm depth | 800 |
| PVO - 15 | 38 cms dia x 45 cm depth | 1000 |
| PVO - 18 | 45 cms dia x 60 cm depth | 1200 |
| PVO - 20 | 50 cms dia x 60 cm depth | 1250 |
| PVO - 24 | 60 cms dia x 70 cm depth | 1500 |