

WATER STILL (WALL MOUNTING - Manesty Type)

Fully made of Stainless steel, electrically operated to obtain pyrogen free distilled water, Models also available with detachable condenser tube for condenser tube cleaning. Suitable to work on 220 / 440 V AC supply.

Cat. No.	Approx: Capacity/hr	watts
PWS – 2	2 ltr	1500
PWS – 4	4 ltr	3000
PWS – 6	6 ltr	4000
PWS – 8	8 ltr.	6000
PWS – 10	10 ltr	8000



WATER STILL (TABLE TOP Branstead type)

Heavy duty model with SS 304 quality material Argon welded.

Cat No.	Approx. Capacity /hr.	Watts
PWS - 2 TT	2 ltr.	1500
PWS - 4 TT	4 ltr.	3000
PWS - 6 TT	6 ltr.	4000
PWS - 8 TT	8 ltr.	6000
PWS - 10 TT	10 ltr.	8000



Automatic water level cut off system can be supplied at an extra cost.