

X-RAY FILM DRIER

This is an equipment for easy and fast drying of X-ray film especially in the monsoon season, equipped with thermostatically controlled heaters and fan for forced air circulation. Provided with hanger slots.

| Cat.No. | Specifications | No.of hanger slot. |
|---------|----------------|--------------------|
| PXD-1 | Single walled | 16 nos. |
| PXD-2 | Double walled | 16 nos. |
| PXD-3 | Single walled | 24 nos. |
| PXD-4 | Double walled | 24 nos. |

